



Technical bulletin

2019-1

This information goes out to medical technicians and resellers of Decapus III.

BATTERY INFORMATION TO END USERS

The battery in Decapus III as a rechargeable Li-ion battery 11.1 V, capacity 2.6 Ah.

When Decapus® III is used for the first time after delivery, the battery must be **fully** charged, which is done by connecting the system to the mains. A flashing green light on the control panel indicates that the battery take charge. This partly needs to be done because the battery on delivery is set in a sleep mode and must be activated. The same applies when changing the battery.

This information can be found in the updated service manual (p.16).

STORAGE CONDITIONS

If the system does not go directly to an end user, but is stocked, attention must be paid to the battery.

After about 4-5 months from production date, the battery must be charged for conservation, which is done by connecting the Decapus III system to the mains. The status of the battery charge level is displayed on the patient module's control panel (4 levels). Charge the battery to level 3, thereafter the system can be stored again.

The same applies to batteries as spare parts and we therefore do not recommend that your stock levels are to high for this part.



Johanna Genberg
Managing Director



May 15th 2019 | TB125_V1

Training certificate

Training	Technical training
Product	EC Sense, EC View, EC Store
Type of training	Functionality, configuration, troubleshooting
Catalogue no.	CN TR-2-1, CN TR-2-2
Trainer	Mattias Carlsson
Training duration	6 h
Date	2019-05-20
Customer	Södersjukhuset

Wajih Habrah

Cardiolex Medical AB hereby ensures that the person above has undergone training in order to perform service, installation, configuration and troubleshooting on above listed product/s.



Mattias Carlsson

Training certificate

Training	Technical training
Product	Decapus III
Type of training	Installation, functionality, troubleshooting
Catalogue no.	QN TR-2
Trainer	Mattias Carlsson
Training duration	6 h
Date	2019-05-21
Customer	Södersjukhuset

Wajih Habrah

Cardiolex Medical AB hereby ensures that the person above has undergone training in order to perform service, installation and troubleshooting on above listed product/s.



Mattias Carlsson